

Glucagon ELISA

Measuring Glucagon in Human Plasma



This cutting-edge kit is designed to accurately measure glucagon levels in human plasma using our advanced ELISA (Enzyme-Linked Immunosorbent Assay) technology.



High Sensitivity



Short Incubation Time



Low Cross-Reactivity




High Precision

The Glucagon ELISA kit is **easy to use** (room temperature incubation & no shaking required), **patient friendly** (no biotin is used), and **end-user friendly** (samples collected in regular EDTA tubes without the addition of inhibitors).

Sensitivity

IBL-America's kit ensures the lowest **LoD**, providing the **highest sensitivity** rating among the kits tested.



LoB: 96 replicates of each lot of Cal A (4 lots total) across 6 tests (3 kit lots, 2 tests for each kit lot) were run.


LoD: 96 replicates of each low-value sample (4 samples in total) across 6 tests (3 kit lots, 2 tests of each kit lot) were run.

LoQ: Lowest glucagon concentration that produces a sample CV% of less than 20%.

Parameter	IBL-America	A	B
LoB	0.80 pg/mL	No Data	No Data
LoQ	5.33 pg/mL	No Data	5.25 pg/mL

Incubation


Time



	IBL-America	A	B
Shaker	No	Yes	Yes
Temperature	Room	Room	4° Overnight

Cross-Reactivity

No detectable cross-reactivity for 21 cross-reactants (including oxyntomodulin and glicentin)



For comprehensive product details, visit: IBL-America.com.



Revision 1.1 – Dec 20, 2024

Glucagon ELISA


Measuring Glucagon in Human Plasma

SINCE 1997




Repeatability

n = 80 results per sample (20 testing days x 2 runs x 2 replicates.) CV ≤ 5.6%



Reproducibility

n = 75 results per sample (3 conditions x 5 testing days x 5 replicates.) CV ≤ 5.6%



Lot-to-Lot Precision

n = 75 results per sample (3 lots x 5 testing days x 5 replicates.) CV: 2.5 - 11.8%


IBL-America	A	B
36 potential interferences were tested.	No interference for hemoglobin, Biotin, intralipids, and bilirubin at two times their physiological concentration.	Interference 36% for oxyntomodulin 1050 pg/mL, hemoglobin > 50 mg/dL, interfere test.
No interference by EP07/EP37 and published literature		

Reference Range

IBL-America	A	B
6.7 - 66.6 pg/mL, n = 124, fasting	13 - 159 pg/mL, undifferentiated population (diabetic and non-diabetic), n = 1533	≤5.2 - 62.7 pg/mL, n=121, fasting

Commutability

Correlation between IBL's Glucagon ELISA and National Reference Laboratory LUMINEX: $r^2 = 0.94$ ($r = 0.97$), slope 0.94 shows **significant**



Kit Accelerated Stability

IBL-America	A	B
At 37°C, stable for up to 11 days	No Data	No Data

In-Use Stability

IBL-America	A	B
Up to 3 months	No Data	No Data

For comprehensive product details, visit:
IBL-America.com.

